



INSTRUCTIONS FOR TESTING, CONSULTATIONS, AND SUBMITTING DIAGNOSTIC SAMPLES

√ Test	Test	Instructions	Cost US\$
	THYROID ANTIBODY PROFILE (D8T) Test includes; T3, T4, free T3, free T4, T3 autoantibodies and T4 autoantibodies.	If the patient is already on thyroid medication the blood should be drawn 4-6 hours after giving thyroid medication. Collect 3 to 5 mL of blood. Add to redtop tube or serum separator tube (SST). Allow to clot and centrifuge to separate serum. DO NOT SEND UNSEPARATED BLOOD. Transfer serum to empty redtop or plastic tube. 1.0 mL serum is REQUIRED for the Thyroid Antibody Profile. Test includes T3, T4, free T3, free T4, T3 autoantibodies and T4 autoantibodies.	\$65.00
	THYROID ANTIBODY PROFILE + THYROGLOBULIN AUTOANTIBODY (D8T + TgAA)	Same as above but send extra 0.5 mL separated serum	\$85.00
	Thyroid Ab Profile, TgAA and OFA REGISTRY PANEL	See above for guidelines except – 2.0 mL serum is REQUIRED to complete all testing. For OFA, include completed & signed OFA form, and submit separate check for \$15.00 made out to OFA.	\$ 85.00 + \$ 15.00
	PROFILE 7200 includes CBC, Differential & Chemistries & D8T	Send 1 mL whole blood in lavender-top tube (EDTA) AND 1.5 mL separated serum as described above for D8T Thyroid Antibody Profile. Include ice pack when shipping, but do not freeze sample. Without the CBC, deduct \$6.	\$95.00
	TgAA (thyroglobulin autoantibody), Add -On	<u>Note:</u> Good for genetic screening. Send extra 0.5 mL separated serum. Allow 3 - 5 days testing time.	\$22.00
	T4 and freeT4 <u>only</u> (44)	Send 0.5 mL separated serum	\$54.00
	CBC, Differential & Chemistries (D2)	Send 1 mL whole blood in lavender-top tube (EDTA) AND 1.5 mL separated serum as described above for D8T Thyroid Antibody Profile. Include ice pack for shipping, but do not freeze. Minus CBC, deduct \$6.	\$55.00
	DISTEMPER & PARVO TITERS	Send 1.0 mL separated serum.	\$36.50
	RABIES VACCINE TITER	Send 1.5 mL serum in red top or plastic tube. <u>Note:</u> Allow 14 - 21 days testing time.	\$75.00
	VON WILLEBRAND TEST wWD	Take 1 mL separated plasma from a citrate tube (light blue) and transfer plasma to empty redtop or plastic tube. Include ice pack when shipping, but do not freeze sample. Clearly label the tube as "PLASMA". <u>Note:</u> Allow 4-7 days testing time.	\$57.00
	PHENOBARBITAL PROFILE	Blood sample should be drawn 4-6 hours after giving medication. 1.0 mL serum is required. <u>Please use red top tube, NOT an SST tube.</u>	\$47.00
	BROMIDE PROFILE	Blood sample should be drawn 4-6 hours after giving medication. 1.0 mL serum is required. <u>Please use red top tube, NOT an SST tube.</u>	\$56.00
	SERUM BILE ACIDS PRE & POST	1st tube is fasting, 2 nd tube is 2-4 hours post-fatty meal. 1.0 mL separated serum in each tube. Label tubes clearly 1 st & 2 nd.	\$47.00
	URINE BILE ACID	Freshly collected urine, 3 mL.	\$45.00
	PANLEUKOPENIA VACCINE TITER – For Felines	Send 1.0 mL separated serum.	\$36.50
	CONSULTATION REQUEST only	(Attach all lab and/or medical information)	\$30.00

- **LABEL TUBES CLEARLY WITH PETS NAME AND OWNERS LAST NAME**
- Ship samples in a **sturdy cardboard box** with plenty of padding.
- **Do not** freeze samples. Blue & lavender top tubes require an ice pack with an overnight courier.
- Send by US Priority Mail, US Express Mail, Fed-Ex, UPS or other courier.
- Samples received in our office Monday – Saturday.
- **HEMOPET** will process and submit samples to Antech Diagnostics for analysis. **Dr. Dodds** will follow with interpretations and recommendations.
- Send check, money order or credit card info along with samples to:

HEMOPET/HEMOLIFE
11330 Markon Drive
Garden Grove, CA 92841

Any questions, please call our office – 714-891-2022 Monday to Friday, 8:30 AM – 5:30 PM.
Download additional forms from <http://www.hemopet.com/services.html>